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- for animals -

Content:

One sterile foldable intraocular lens made of hydrophilic acrylate (25% - 28% water content) with incorporated UV absorber and hydrophobic surface.
See label on the cardboard box for type and diopter.

Description:

Each Acrivet intraocular lens is a precision-made optical device finished by the finest workmanship. Manufacture of the lenses complies with international standards. Acrivet acrylic intraocular lenses consist of ultrapure, biocompatible acrylate with chemically bound UV absorber.
The haptics of the 3-piece acrylic IOLs are made of Loop, a blue, steam-sterilizable monofilament.
The tolerance of refractive power of an Acrivet intraocular lens is ± 0.25 D for a dioptric range of < 25.0 D and ± 0.5 D for a dioptric range > 25.0 D and ≤ 30.0 D. IOLs with a refractive power of > 30.0 D have a tolerance of ± 1.0 D.

Sterilization:

Acrivet acrylic intraocular lenses are terminally sterilized by autoclaving according to a validated procedure in their original vials.

Packaging:

Each intraocular lens is packaged sterile in a lens holder in a vial filled with ultrapure water, which is packaged in a sterile peel pouch.

Do not use the intraocular lens if the package is damaged.

Handling:

Examine the label on the lens package for IOL type, dioptic power and expiration date. To remove the IOL, open the Medipeel pouch at the defined spot and remove the lens vial in a clean environment. Verify the dioptic power of the lens again. Open the lens vial and remove the lens ensuring a sterile condition. Rinse the lens thoroughly using sterile intraocular irrigating solution such as BSS (Balanced Salt Solution).

Avoid desiccation of the lens – especially between time of removal from the vial and time of implantation – in order to maintain foldability.

Acrivet acrylic intraocular lenses are foldable in hydrated condition. They may be implanted either using a suitable injector or a folding forceps.

Indication:

- Aphakia following surgical removal of the crystalline lens in the animal's eye

Contraindications:

1. Absolute contraindications:

Predominantly chronic cases with potentially unfavourable results according to experience.

- Chronic or recurring uveitis
- Severe corneal dystrophy
- Uncontrolled glaucoma
- Acute eye disease or infection

2. Relative contraindications

Clinical cases which may deteriorate due to IOL implantation and cases with an increased risk caused by implantation. The individual evaluation of each single case is left to the surgeon.

3. Surgical contraindications

- Hyphema
- Excessive vitreous loss
- Zonular damage
- Posterior capsular rupture
- Presence of, or predisposition to, retinal detachment

Side effects / interactions:

No side effects or drug interactions have been reported.

Precautions:

1. Do not resterilize intraocular lenses.
2. Do not implant intraocular lenses from damaged Medipeel pouches.
3. Never use intraocular lenses beyond the expiration date.
4. Store intraocular lenses away from moisture and light, between 2°C (35.6°F) and 50°C (77°F). Do not use any lenses from moist packages.
5. Do not use the storage medium as intraocular irrigation solution.
6. If the Intraocular lenses will be stored at lower temperatures ($< 3^{\circ}\text{C}$, e.g. while transport) reversible short time opacification can occur. Intraocular lenses should have room temperature at the time of implantation in order to avoid this opacification.
7. In case of implantation with a cartridge, the lens may show expression furrows. These furrows are reversible and therefore give no reason for explantation.
8. The surgeon should receive training for this type of surgery prior to attempting implantation of this lens.

Return of unsterile products or damaged package:

Only returns including the original packaging which identifies the lens and its serial number will be accepted. Request a Return Merchandise Authorization No. (RMA-No.) in advance (Customer-Service@skin-vision.eu).

Contaminated products should not be sent back.

Disclaimer of liability:

Valeant Med Sp. z o.o. shall not be liable for any misuse of the product or for the selection and use of the lens for an individual patient or the condition of the patient.

Used symbols:

| | | | |
|---|---------------------|---|------------------------------------|
|  | Steam-sterilized |  | Store away from sunlight |
|  | Serial number |  | Expiry date |
|  | Do not reuse |  | Observe instructions for use |
|  | Store in dry places |  | Storage temperature |
|  | Do not re-sterilize |  | Do not use when package is damaged |

BAUSCH + LOMB

 Valeant Med Sp. z o.o.
ul. Ryzowa 31
02-495 Warszawa, Poland
Tel: +48 22 578 16 00
Fax: +48 22 578 16 18

Acrivet is a trademark of Valeant Pharmaceuticals International, Inc. or its affiliates.

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©Bausch & Lomb Incorporated.